

# **EXHIBIT B20**

## **Part 3**

RayBio # ELH-CA125

Test unconcentrated media vs concentrated using Amecim  
Ultra-15 filter MW cutoff 10000

<u>Weigh tubes</u>	<u>Volume empty</u>	<u>empty sample reservoir w/ media</u>	<u>empty centrifuge tube</u>	<u>media + reservoir</u>
266	8.36 mL	25.0 g	10.99 g	11.59 g
338	8.853 mL	24.9 g	10.98 g	11.53 g
				19.17 g
				19.7 g

Spin tubes for 25 min at 4000xg

move the retentate by pipetting into new container

<u>weigh retentate</u>	<u>tube weight</u>	<u>tube + retentate</u>	<u>filtrate</u>
266	1.0025 g	1.3514 g	19.29 g
338	1.0088 g	1.7104 g	19.4797 g


$$\% \text{ retentate recovery} = 100 \times \frac{W_r \times C_r}{W_o \times C_o}$$

$$\% \text{ filtrate recovery} = 100 \times \frac{W_f \times C_f}{W_o \times C_o}$$

$$\% \text{ recovery} = \% \text{ retentate recovery} + \% \text{ filtrate recovery}$$



1. All reagents & samples to room temperature
2. Assay diluent (Item E2) should be diluted 5x w/ dd H<sub>2</sub>O
  - Stable 1 mo. at 4°C
3. Prep. of standard
  - Spin vial C
  - Add 400  $\mu$ L 1x Assay Diluent into vial C = 1000 u/mL
  - mix gently



4. Prepare Wash buffer by diluting 20x
  - Stable 1 mo. at 4°C

5. Spin Hem F, detection antibody
  - Add 100  $\mu$ L of 1x Assay Diluent
  - Stable 5 days at 4°C
  - Dilute it ~~80x~~ <sup>80x</sup> and will be used in assay

6. Spin HRP - strep. Vial (Item G) & mix
  - Add 15  $\mu$ L to tube w/ 12 mL of assay diluent (800x)
  - do not save!

Go to pg 10



### Assay Procedure

1. All samples to RT.
2. Label 8 strip/wells
3. Add 100ul standard / samples
  - Incubate 2.5 hr, gentle shaking, RT
4. Discard solution, wash 4x
  - decant, blot after washes
5. Add 100ul of 1x Antibody mix
  - Incubate 1 hr, RT, shaking
6. Discard & wash as in step 4 (4x)
7. Add 100ul of Streptavidin solution
  - Incubate 45 min, RT, shaking
8. Discard & wash 4x
9. Add 100ul TMB substrate (Item H)
  - Incubate 30 min, RT, shaking
10. Add 50ul STOP solution (Item I) to each well
  - Read at 450 nm



	1	2	3	4	5	6	7	8	9	10	11	12
A	Stand	A		Stand	I							
B		B		Blank								
C		C		266 orig								
D		D		338 orig								
E		E		266 conc.								
F		F		338 conc.								
G		G										
H		H										

Standard	Raw1	Raw2	Raw3	Average	Conc (ug/ml)
A	0.5417	0.4288	0.4595	0.44415	2000
B	0.3952	0.3971	0.3944	0.392233333	1500
C	0.3347	0.3186	0.3157	0.323	1000
D	0.2741	0.2863	0.2583	0.266233333	750
E	0.2216	0.214	0.2153	0.216966667	500
F	0.1689	0.1608	0.1584	0.162033333	250
G	0.1374	0.1405	0.1339	0.137233333	125
H	0.1105	0.1072	0.1072	0.1083	25
I	0.1171	0.1195	0.1258	0.1208	5
Blank	1.0783	1.1058	1.0779	1.09185	0

Y = 0.0002x + 0.1222  
 R<sup>2</sup> = 0.98251

	Raw1	Raw2	Raw3	Corr 1	Corr2	Corr3	ug/ml 1	ug/ml 2	ug/ml 3	ug/ml	ug/ml	average	ul to use	diluent to u	ul to use	diluent to use
266 original	0.8548	0.8586	0.8643				3683	3682	3710.5	3683	3682	3.7	95.0	5.0	332.4	17.6
338 original	0.952	0.9891	1.0074				4149	4304.5	4428	4149	4335	4.3	81.3	18.7	284.7	65.3
266 concentrated	4	4	4				19389	19389	19389	19389	19389	19.4	100.0000	0	350.0	0.0
338 concentrated	4	4	4				19389	19389	19389	19389	19389	19.4	100.0000	0.0	350.0	0.0




CA125 ELISA - test levels in media

	1	2	3	4	5	6	7	8	9	10	11	12
A	Stand	1		266 orig								
B		2		338 orig.								
C		3		266 conc								
D		4		338 conc								
E		5		media blank								
F		6										
G		7										
H	Blank											

1/11/18 Test media amounts for CA125 ELISA											
Standard (U/ml)	OD1 (450 nm)	OD2	Corrected OD1	Corrected OD2	Average						
400	2.3856	2.3921	2.31895	2.32545	2.3222						
133.3	1.1825	1.1458	1.11585	1.07915	1.0975						
44.45	0.3643	0.358	0.29765	0.29135	0.2945						
14.81	0.1593	0.1562	0.09265	0.08955	0.0911						
4.94	0.1049	0.1009	0.03825	0.03425	0.03625						
1.65	0.082	0.0861	0.01535	0.01945	0.0174						
0.55	0.0769	0.0776	0.01025	0.01095	0.0106						
Blank	0.0661	0.0672	0.06665								



$y = 0.0065x + 0.004$   
 $R^2 = 0.9982$

Samples	OD1	OD2	OD3	Corrected OD1	Corrected OD2	Corrected OD3	Corrected for Media OD1	Corrected for Media OD2	Corrected for Media OD3	CA125 U/ml	CA125 U/ml	CA125 U/ml	Average	Standard Deviation
Original 266	0.096	0.0925	0.0939	0.02935	0.02585	0.02725	0.0155	0.012	0.0134	1.769230769	1.230769231	1.446153846	1.482051282	0.271019698
Original 338	0.1179	0.1155	0.1172	0.05125	0.04885	0.05055	0.0374	0.035	0.0367	5.138461538	4.769230769	5.030769231	4.979487179	0.189882139
Concentrated 266	0.0849	0.0799	0.0843	0.01825	0.01325	0.01765	0.0044	-0.0006	0.0038	0.061538462	-0.70769231	-0.030769231	-0.22564103	0.42001221
Concentrated 338	0.2216	0.2312	0.2172	0.15495	0.16455	0.15055	0.1411	0.1507	0.1367	21.09230769	22.56923077	20.41538462	21.35897436	1.101406667
Plain Media	0.0799	0.0885	0.0811	0.01325	0.02185	0.01445	0.01385							

266 = T0V120 unit

338 = Normal Ov. Epi cells

Proceed using unconcentrated media




# Protein levels for CA125 assay

Re-did standard and original media - media was too concentrated

• diluted media by 50%, remeasured

• Also used 25 ul of the standard and samples




1/12/2018

Samples	OD1	OD2	OD3	Corrected OD1	Corrected OD2	Corrected OD3	CA125 U/ml	CA125 U/ml	CA125 U/ml	Average	Standard Deviation
TOV112D 72 hr control	279	0.104	0.0885	0.0969	0.03735	0.02185	0.03025	5.13076923	2.74615385	4.03846154	4.58461538
TOV112D 72 hr 1000 ug/ml Talc	281	0.0806	0.0742	0.0838	0.01395	0.00755	0.01715	1.53076923	0.54615385	2.02307692	0.77237818
OV90 72 hr control	297	0.1059	0.1046	0.1321	0.03925	0.03795	0.06545	5.42307692	5.22307692	9.45384615	1.77692308
OV90 72 hr 1000 ug/ml Talc	299	0.1062	0.0906	0.1316	0.03955	0.02395	0.06495	5.46923077	3.06923077	9.37692308	0.34811411
TOV-21G 72 hr control	333	0.0724	0.0716	0.0942	0.00575	0.00495	0.02755	0.26923077	0.14615385	3.62307692	0.14142136
TOV-21G 72 hr 1000 ug/ml Talc	335	0.0761	0.077	0.0856	0.00945	0.01035	0.01895	0.83846154	0.97692308	2.3	1.69705627
Normal Ov Epithelial 72 hr control	351	0.1003	0.0943	0.1001	0.03365	0.02765	0.03345	4.56153846	4.53076923	4.54615385	0.08702853
Normal Ov Epithelial 72 hr 1000 ug/ml talc	353	0.1106	0.092	0.1331	0.04395	0.02535	0.06645	6.14615385	3.28461538	9.60769231	0.09790709
b		0.0661	0.0672		0.06665					4.71538462	0.02175713
										2.02341325	

The other proteins in media may be interfering. Try lysate.



Standard ID	Concentration n (ug/ml)	OD1	OD2	OD3	Average	Corrected Avg
A	2000	0.5869	0.5848	0.5562	0.57596667	0.4874
B	1500	0.4457	0.4211	0.3533	0.4334	0.34483333
C	1000	0.3213	0.2774	0.2593	0.286	0.19743333
D	750	0.2473	0.2199	0.2128	0.22666667	0.1381
E	500	0.181	0.1759	0.1713	0.17606667	0.0875
F	250	0.1328	0.1252	0.1256	0.12766667	0.0393
G	125	0.1105	0.112	0.1116	0.11136667	0.0228
H	25	0.1004	0.0941	0.0953	0.0966	0.00803333
I	5	0.0845	0.0846	0.0941	0.08773333	0
J (BLANK)	0	0.0881	0.0884	0.0892	0.08856667	0



											using 100			
	ID	OD1	OD2	OD3	Corr OD1	Corr OD2	Corr OD3	ug/ul 1	ug/ul 2	ug/ul 3	Average	ug per well	x 3.5 wells	diluent
TOV112D 72 hr control	279	0.2429	0.2301	0.2326	0.1356	0.1228	0.1253	7.605	6.965	7.09	7.0275	14.23	49.8	300.2
TOV112D 72 hr 1000 ug/ml Talc	281	0.1892	0.1923	0.1878	0.0819	0.085	0.0805	4.92	5.075	4.85	4.94833333	20.21	70.7	279.3
OV90 72 hr control	297	0.3198	0.331	0.3306	0.2125	0.2237	0.2233	11.45	12.01	11.99	12	8.33	29.2	320.8
OV90 72 hr 1000 ug/ml Talc	299	0.2689	0.2759	0.2784	0.1616	0.1688	0.1711	8.905	9.255	9.38	9.3175	10.73	37.6	312.4
TOV-21G 72 hr control	333	0.2991	0.2968	0.2948	0.1918	0.1895	0.1875	10.415	10.3	10.2	10.305	9.70	34.0	316.0
TOV-21G 72 hr 1000 ug/ml Talc	335	0.1658	0.1616	0.1634	0.0585	0.0543	0.0561	3.75	3.54	3.63	3.64	27.47	96.2	253.8
Normal Ov Epithelial 72 hr control	351	0.2359	0.23	0.23	0.1286	0.1227	0.1227	7.255	6.96	6.96	6.96	14.37	50.3	299.7
Normal Ov Epithelial 72 hr 1000 ug/ml talc b	353	0.2432	0.2474	0.2344	0.1359	0.1401	0.1271	7.62	7.83	7.18	7.725	12.94	45.3	304.7
		0.1097	0.1049		0.1073									
											using 100			
	ID	OD1	OD2	OD3	Corr OD1	Corr OD2	Corr OD3	ug/ul 1	ug/ul 2	ug/ul 3	Average	ug per well	x 3.5 wells	diluent
TOV112D 72 hr unt	278	0.6686	0.6925	0.653	0.56486667	0.58876667	0.54926667	29.06833333	30.26333333	28.28833333	29.27583333	3.42	12.0	338.0
OV90 72 hr unt 72 hr	296	0.34	0.3407	0.3366	0.23626667	0.23696667	0.23286667	12.63833333	12.67333333	12.46833333	12.59333333	7.94	27.8	322.2
TOV-121G unt	332	0.4392	0.4456	0.4413	0.33546667	0.34186667	0.33756667	17.59833333	17.91833333	17.70333333	17.81083333	5.61	19.7	330.3
Normal Ov Epithelial 72 hr unt	350	0.2281	0.2347	0.2316	0.12436667	0.13096667	0.12786667	7.04333333	7.37333333	7.21833333	7.29583333	13.71	48.0	302.0
b		0.1032	0.104	0.104	0.10373333									

	1	2	3	4	5	6	7	8	9	10	11	12
A											279	
B											281	
C											297	
D											299	
E											333	
F	BLANK								BLANK		335	
G	279								BLANK		351	
H	281								BLANK		353	

	1	2	3	4	5	6	7	8	9	10
A	Stand	A			I					
B	B				PBS Blank					
C	C				278					
D	D				296					
E	E				332					
F	F				350					
G	G				Lysis Buffer Blank					
H	H									



	1	2	3	4	5	6	7	8	9	10	11	12
A	279											
B	281											
C	297											
D	299											
E	333											
F	335											
G	351											
H	353											

	1	2	3	4	5	6	7	8	9	10	11	12
A											278	
B											278	
C											296	
D											296	
E											332	
F											332	B
G											350	332
H											350	350

1/17/2018	Lysate	OD1	OD2	OD3	Corrected OD1	Corrected OD2	Corrected OD3	CA125 U/ml	CA125 U/ml	CA125 U/ml	Average	Standard Deviation
TOV112D 72 hr control	279	0.1637	0.1749	0.171	0.0807	0.0919	0.088	11.8	13.5230769	12.9230769	13.2230769	0.42426407
TOV112D 72 hr 1000 ug/ml Talc	281	0.1188	0.1175	0.1221	0.0358	0.0345	0.0391	4.89230769	4.69230769	5.4	4.99487179	0.36482413
Ov90 72 hr control	297	0.146	0.1479	0.1689	0.063	0.0649	0.0859	9.07692308	9.36923077	12.6	9.22307692	0.20669275
Ov90 72 hr 1000 ug/ml Talc	299	0.1371	0.13	0.1506	0.0541	0.047	0.0676	7.70769231	6.61538462	9.78461538	6.29230769	0.77237818
TOV-21G 72 hr control	333	0.1374	0.1269	0.1289	0.0544	0.0439	0.0459	7.75384615	6.13846154	6.44615385	6.43076923	0.19581419
TOV-21G 72 hr 1000 ug/ml Talc	335	0.1	0.0889	0.0907	0.017	0.0059	0.0077	2	0.29230769	0.56923077	9.76153846	0.29372128
Normal Ov Epithelial 72 hr control	351	0.1491	0.1518	0.1547	0.0661	0.0688	0.0717	9.55384615	9.96923077	10.4153846	9.83076923	0.91379953
Normal Ov Epithelial 72 hr 1000 ug/ml talc	353	0.1551	0.1362	0.1467	0.0721	0.0532	0.0637	10.4769231	7.56923077	9.18461538	9.83076923	0.91379953
b		0.0868	0.0876	0.0746	0.083							

1/18/2018	Lysate	OD1	OD2	OD3	Corrected OD1	Corrected OD2	Corrected OD3	CA125 U/ml	CA125 U/ml	CA125 U/ml	Average	Standard Deviation
TOV112D 72 hr unt	278	0.1998	0.1904	0	0.0476	0.0382	0.0468	6.70769231	5.26153846	6.4	5.98461538	1.02258519
Ov90 72 hr unt 72 hr	296	0.1978	0.1971	0.199	0.0456	0.0449	0.0468	4.81538462	3.50769231	6.58461538	4.16153846	0.07614996
TOV-121G unt	332	0.1875	0.179	0.2328	0.0353	0.0268	0.0806	14.9384615	20.1846154	11.7846154	13.3615385	2.230106
Normal Ov Epithelial 72 hr unt	350	0.2533	0.2874	0.2328	0.1011	0.1352	0.0806	14.9384615	20.1846154	11.7846154	13.3615385	2.230106
b		0.1522										

Results: Lysate protein measurements may be affected by talc.  
Repeat protein measurements, have control w/ talc in it.